



To determine the oil or finish content of wool or synthetic samples by solvent percolation and evaporation. Suitable for fiber, yarn or fabric. Available with 1, 2 or 3 extraction units. To be fully compliant with Edana method ERT 155, it needs a large extraction tube- 12 mm (OD), 10 mm (ID). (Full title of new method is: DETERMINATION OF ACETONE EXTRACTABLE FINISH ON NONWOVENS).

SIZE

WEIGHT 4.4 kg / 9.7 lbs

DIMENSIONS 170 mm (w) x 170 mm (d) x 380 mm (h) /
7 in (w) x 7 in (d) x 15 in (h)

ORDERING INFORMATION

105265 Rapid Oil Extraction Apparatus Single

105266 Rapid Oil Extraction Apparatus Double

105267 Rapid Oil Extraction Apparatus Triple

105268 Foil Residue Trays (Pack of 50)